
TERMAKS

SERIES TS8000 Drying Ovens / Sterilizers

USERS MANUAL / *BRUKSANVING*

ENGLISH / <i>NORSK</i>	Page	Side
General Description / <i>Generell Beskrivelse</i>	2	9
Installation / <i>Installasjon</i>	2	9
Flowchart, Basic operations / <i>Flytskjema, Hovedfunksjoner</i>	3	10
Temperature / <i>Temperatur</i>	4	11
Check The Temperature Set / <i>Vis innstilt temperatur</i>	4	11
Change The Temperature Set / <i>Endre innstillingsverdi</i>	4	11
Safety Thermostat / <i>Sikkerhetstermostat</i>	4	11
Datalogging / <i>Datalogging</i>	4	11
Alarm / <i>Alarm</i>	5	12
Change Alarm Limit / <i>Endre Alarmgrense</i>	5	12
Time / <i>Tid</i>	5	12
Check The Time Set / <i>Vis Innstilt Tid</i>	5	12
Change The Time Set / <i>Endre Tidsinnstillingen</i>	5	12
On/Off / <i>Av/På</i>	6	13
Cycle Completed / <i>Syklus Fullført</i>	6	13
Cycle Repeat / <i>Repetert Syklus</i>	6	13
Power Failure/Timer / <i>Strømprudd/Tidsur</i>	6	13
Fan Control / <i>Viftestyring</i>	6	13
Change The Fan Speed / <i>Endre Viftehastighet</i>	6	13
Stop At Door Opening / <i>Viftestopp Ved Døråpning</i>	6	13
Flowchart, Alarm, Calib., Options / <i>Flytskjema, Alarm, Kalib., opsj.</i>	7	14
Calibration / <i>Kalibrering</i>	7	14
Change The Calib. Constant / <i>Endre Kalib. konstant</i>	7	14
Options / <i>Opsjoner</i>	8	15
Program Release / <i>Programversjon</i>	8	15
Troubleshooting / <i>Feilsøking</i>	8	15

GENERAL DESCRIPTION

All control functions are built into a microprocessor based electronic controller.

The main functions are:

- * Main temperature controller
- * Acoustic and visual alarm system
- * Timer, counts only when the temperature is at setpoint. Acoustic warning at time-out.
- * Fan speed control and fan stop at door opening
- * Independent safety thermostat for material protection
- * The cabinet has fused power inlet and a thermal temperature switch for maximum safety.

Options available for the user.

Hours instead of minutes, expand to 250°C working range, delayed start, and disabling of the acoustic alarm warning.

INSTALLATION

All three sizes are table models and are stackable. The two smallest models can be hung on a wall. Separate brackets for this purpose is delivered with those models. The biggest model can be placed on an optional floor stand.

Avoid flammable materials behind the cabinet, such as wall paper. Also avoid storage of flammable materials on the top of the cabinet.

The cabinet should not be placed in a room with high humidity or corrosive atmosphere


The cabinet shall be connected to a grounded 230 VAC line.

Turn on the cabinet with the power switch.

If the timer is not to be used, enter 000 as time set to prevent accidental activating of the timer.

If necessary, enter a new option code

Flowchart Basic operations Temperature - Time - Fan Speed

TEMPERATURE

Normally, the °C symbol is lit, indicating that the system is in the temperature display mode. The system returns automatically to this mode within a few minutes if the cabinet is left in an other display mode, and this is the default mode when turning on the cabinet. One or two hits on the *Select* button returns the system to the temperature display mode from the fan speed or time mode.

Check the temperature set

The temperature set is displayed simply by pressing the *Set* button

Change the temperature set

Press the *Set* button for 3 seconds, then the °C symbol will start flashing.

Release the *Set* button and enter the new set value with the *Up* or *Down* buttons.

The set value will change step by step at short hits. If they are continuously pressed, the value will first change in a slow rate, then fast.

Maximum set value is 200°C, but can be expanded to 250°C. (Read the option chapter.)

When the new set is entered, it can be stored in the memory by hitting the *Set* button, or it will automatically be stored after a 8 second delay. The new set value is stored in the memory when the °C symbol stops flashing.

Safety thermostat


It is very important to adjust the safety thermostat to the same value as entered in the display. The trip temperature for this thermostat is normally about 10-15°C above the temperature indicated on the scale. In some cases, the materials inside the cabinets is very sensitive to overtemperature. In such cases, the safety thermostat should be carefully adjusted down to the trip point after temperature stabilisation. Then turned a little bit up again. This will give maximum protection of such sensitive materials. If this thermostat is adjusted too close to the working temperature, it may affect the control accuracy and the timer function.

Datalogging

The maximum, minimum and the average temperature values during a cycle are stored in the controller.

These values are displayed simply by pressing the appropriate buttons. This can be an useful aid in the evaluation of the temperature cycle.

Do not press the Average button for more than 2 seconds to avoid changes in the timer on/off state or a restart of the cabinet.

The calculations are done on a minute basis, and starts about 32 minutes after the initial temperature stabilising. In this initial period, "- - -" will be displayed if any of those buttons are pressed. Therefore, the datalogging will never start at cycles less than 32 minutes.

If the timer function is turned on, the data collection stops when the cycle is completed.

The values are kept in the memory as long as the green light, the cycle completed indicator, is lit.

Therefore, all those data can be displayed at any time before the cabinet are restarted.

An average value very close to the temperature setpoint indicates that the temperature must have been correct during most of the cycle.

A large deviation of the max or min values in this case, indicates a temporary disturbance, for example a short door opening.

Obviously, in very long cycles, a more serious deviation will still have very little influence on the average value.

The values are cleared in the following cases:

- * Initial start
- * Change of setpoint
- * Restart after a ended time cycle.
- * After a power failure, except if the power fails after a ended time cycle.

Alarm

If the temperature should happen to drift outside the temperature limit, the display will start to flash between actual temperature and three dashes. If the alarm condition lasts for more than about 4 minutes, an acoustic warning will be turned on as well. This warning can be turned off by a hit on the *Set* button. The limit is common for over- and undertemperature.

The alarm is disabled until the temperature has stabilised after a setpoint change or a restart of the cabinet.

Change alarm limit (see also page 7)

New alarm limit can be entered by pressing the *Set* button for 3 seconds until the °C symbol starts flashing. Then release the button and press it in again for 6 seconds until °C and the *Time* symbol is flashing simultaneously. Actual alarm limit is displayed and may be changed by the *Up/Down* buttons. When the new limit is entered, it can be stored in the memory by hitting the *Set* button, or it will automatically be stored after a 8 second delay. The new set value is stored in the memory when the symbols stop flashing.

The limit is factory set to 5°C.

TIME

Hit the *Select* button one or two times until the *Time* symbol is lit. This indicates that the system is in the time display mode, normally minutes.

If the timer is on, the remaining time of the cycle is displayed.

If the timer is off, - - is displayed, indicating no relevant data is available.

Check the time set

The time set is displayed simply by pressing the *Set* button

Change the time set

Press the *Set* button for 3 seconds, then the Time symbol will start flashing.

Release the *Set* button and enter the new set value with the *Up* or *Down* buttons.

The set value will change step by step at short hits. If they are continuously pressed, the value will first change in a slow rate, then fast.

When the new set is entered, it can be stored in the memory by hitting the *Set* button, or it will automatically be stored after a 8 second delay. The new set value is stored in the memory when the Time symbol stops flashing.

This time set means the time from the working temperature is reached to the cabinet is turned off at the end of cycle.

If the time set is changed in *Timer on* mode, the timer will count down from that new time, and future cycles will use this new time set.

If the set value 000 is entered, the timer will automatically turned off, and it is not possible to turn the timer on again as long as 000 is entered. This should always be done in cases where the timer not shall be used in order to prevent accidental activating of the timer.

On/off

The timer is turned on by pressing the *Timer on/off* button for 3 seconds until the yellow timer on indicator is lit.

From the on state, it can be turned off again by pressing the *Timer on/off* button for 3 seconds until the yellow timer on indicator is turned off.

The *Timer on/off* button can be operated in any display mode.

In the off state, yellow and green lights off, the cabinet will work continuously.

Cycle completed

When the timer is on and the temperature has been at set temperature for the actual time set, the green cycle completed light will be lit, and the heater is turned off. A acoustic signal will also be heard for a few seconds. The fan will remain on for about 8 minutes in order to prevent a temperature overshoot.

Cycle repeat

In the Cycle completed state, the timer can be turned on again by pressing the *Timer on/off* button for 3 seconds until the yellow timer on indicator is lit.

The cabinet will then repeat the previous cycle if none of the settings are changed since last cycle.

Power failure/Timer

In the case of a power failure during a cycle, the cabinet will repeat the cycle from start.

A power failure in the cycle completed mode will not affect the cabinet.

All logged data are kept in the memory.

FAN CONTROL

Hit the *Select* button one or two times until the *Fan symbol* is lit. This indicates that the system is in Fan speed display mode. The numbers 0,1,2,3 or 4 will appear in the display.

0 means that the system has stopped the fan for some reason. It is not possible for the user to enter this stop code.

Change the fan speed

Press the *Set* button for 3 seconds, then the *Fan symbol* will start flashing.

Release the *Set* button and enter the new speed code with the *Up* or *Down* buttons.

The set value will change step by step at short hits.

The codes 1,2,3 and 4 can be entered by the user. 1 is the lowest while 4 is full speed.


When the new speed code is entered, it can be stored in the memory by hitting the *Set* button, or it will automatically be stored after a 8 second delay. The new speed code is stored in the memory when the *Fan symbol* stops flashing.

Fan stop at door opening

In order to keep the hot air inside the cabinet, and to do the work situation for the operator more comfortable, the fan will stop for up to 40 seconds after a door opening. It stops when a significant drop in the temperature occur, but will start again if the temperature comes back to the working temperature, or at last 40 seconds later.

This fan stop routine is disabled for 8 minutes from the last door opening.

Flowchart Alarm limit - Zero calibration - Options

CALIBRATION

A calibration constant may be changed in order to bring the temperature measurement equal to the true value according to a control device

Change the calibration constant

Press the *Set* button for 3 seconds until one of the symbols starts flashing.
Release the button and press again for 6 seconds until the two upper symbols start flashing.
Release the button and press again for 6 seconds until the two lower symbols start flashing.
Release the button and read the old calibration constant. Then calculate the new constant.

Example: Displayed temperature:
True temperature:
Old calibration constant:

180°C
178,5°C
85

$$\text{New constant} = 85 + (180 - 178,5) \cdot 4 = 91$$

If the symbols are still flashing, the new constant can be entered immediately with the *Up* and *Down* buttons. If not, this procedure has to be repeated

When the new constant is entered, it can be stored in the memory by hitting the *Set* button, or it will automatically be stored after a 8 second delay. The new constant is stored in the memory when the flashing stops.

OPTIONS

The user can choose one or more of the following options:

Hours (code 1)

The timer will use hours as time base instead of minutes

250°C (code 2)

Expands the working range from 200 to 250°C

Delayed start (code 4)

The standard timer is replaced by a delayed start timer. The entered time is the start delay.

The delay starts to count down when the timer is turned on by pressing the *Timer on/off* button for 3 seconds until the yellow light is turned on. At time-out, this light will be turned off and the cabinet will start in continuous mode. The entered delay is stored for the next delay period.

Disable acoustic alarm (code 8)

The acoustic warning at alarm conditions is disabled.

Change the option code

Press the *Set* button for 3 seconds until one of the symbols starts flashing.

Release the button and press again for 6 seconds until the two upper symbols start flashing.

Release the button and press again for 6 seconds until the two lower symbols start flashing.

Release the button and press again for 6 seconds until the upper and lower symbols start flashing.

Release the button and enter the new option code with the *Up* or *Down* buttons.

Enter the appropriate code or the sum of the options desired.

When the new code is entered, it can be stored in the memory by hitting the *Set* button, or it will automatically be stored after a 8 second delay. The new code is stored in the memory when the flashing stops.

PROGRAM RELEASE

The software in the microprocessor may be changed from time to time of different reasons. If you have any questions regarding the system, please read out the release number before calling Termaks.

The release number will be displayed when the *Up* and *Down* buttons are pressed simultaneously.

TROUBLESHOOTING

Dead system, no lights

Check the fuses in the power inlet.

Disconnect the cabinet from the power line, unscrew the two screws on the front panel, pull out the controller and check the fuse F1 on the board. If the fuse is blown, check if the fan movement should happen to be blocked of some reason before replacement of the fuse.

Unstable temperature control

Check that the fan is working properly

Check the safety thermostat adjustment


Timer counts down too slow

Check the safety thermostat adjustment.

Check if the hour option code is entered

No Heat, even if the heat indicator is on

Check the thermal switch close to the fan motor in the rear of the cabinet. If this switch is open, the reason for the failure should be fixed before resetting the switch with the red button.

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