UiO : University of	Oslo	
Job Safety Analysis (JSA)	Case and doc. no in ePhorte: 2012/11454-30
		Page : 1 of 2
Prepared by:		Date : 24.09.2012
HSE staff, OPA		Edition: 1

Unit: Nikon XT H 225 ST micro computed tomography (micro-CT) Date: __/__/20__

Job task to which this Job Safety Analysis applies: Non-destructive high-resolution 3D imaging of solid objects with a Nikon XT H 225 ST micro-CT

Responsible employer:
Time of performance of the work:
Persons involved in performance of the work:
•

Description of the work: Scanning of solid materials

Training required: Introduction to the operation of the Nikon XT H 225 ST micro-CT (length of the training variable according to the users experience).

PPE and other protective equipment to be used: eventually gloves (during sample preparation).

Analysis:

Hazards	Preventive measures
Ionization radiation	Micro-CT is regularly serviced and safety checked.
	There should be no exposure to X-rays or to the ionization radiation.
	The operator is not exposed to any radiation under normal working mode. Never deviate from normal moor try to override the interlock.
	If you are a regular user of the lab, use your dosimeter.
Electrocution	Do not remove the protective covers from the micro-C instrument. Do not open other components of the equipment. Do not touch cables and wires. Do not let water or other objects in contact with electronic
	Ionization radiation

UiO: University of	Oslo	
Job Safety Analysis (J	SA)	Case and doc. no in ePhorte:
		2012/11454-30
		Page : 2 of 2
Prepared by:		Date : 24.09.2012
HSE staff, OPA		Edition: 1

Sub-tasks	Hazards	Preventive measures
		components of the equipment.
Introducing and/or removing samples in the micro-CT	Personal injury caused by jamming body parts (hands)	Make sure hands and/or other body parts are clear of t chamber when closing chamber door.
Prolonged use	Eye strain	Do regular breaks (10 min/h during operation)
	Tiredness	Do regular breaks (10 min/h during operation)
	Musculoskeletal strain or injury	Do regular breaks (10 min/h during operation)
		Adjust height of the seat, so that the back and neck are upright

Participants in the analysis and signatures

Role	Name	Date/signature
Lab Manager/d	eputy	/
		/
		/
		•
The Job Safet the work)	y Analysis has been rea	${f d}$ (signatures of users who will perform or be affected b
	y Analysis has been rea	${f d}$ (signatures of users who will perform or be affected b ${f Date/signature}$
the work)	y Analysis has been rea	
he work) Name		